

BIL	Quotation Reference	Description	Advertisement Date	Closing Date (Not Later Than 02.00PM)	Quotation Fee	Requesting Department	Focal Person
2	BME/A2K/(1194)/2023/HTD(B8)	<p>ONE-OFF PREVENTIVE MAINTENANCE SERVICES FOR TWELVE ONE (12) UNITS OF VARIOUS BRAND AND MODELS PHARMACUTICAL REFRIGERATORS AT VACCINE CENTERS BERAKAS HEALTH SCREENING, MINISTRY OF HEALTH</p> <p>PLACE OF SUBMISSION: QUOTATION BOX (GROUND FLOOR) MINISTRY OF HEALTH COMMONWEALTH DRIVE BANDAR SERI BEGAWAN, BB 3910 NEGARA BRUNEI DARUSSALAM</p>	02/12/2023	23/12/2023	\$5.00	JABATAN TEKNOLOGI PENJAGAAN KESIHATAN, KEMENTERIAN KESIHATAN.	<p>HAJI AMIR HISMA BIN HAJI MASRI</p> <p>JURUTERA PERUBATAN HAYAT</p> <p>TEL: 2340228 ext. 511 (HC)</p>

LAMPIRAN - A1

RUJ: BME/A2K/(1194)/2023/HTD (B8)

28/11/2023

ONE-OFF PREVENTIVE MAINTENANCE FOR TWELVE (12) UNITS OF VARIOUS BRANDS AND MODELS PHARMACEUTICAL REFRIGERATORS AT VACCINE CENTERS BERAKAS HEALTH SCREENING, MINISTRY OF HEALTH.

SCOPE OF WORK

ONE-OFF PREVENTIVE MAINTENANCE FOR TWELVE (12) UNITS OF VARIOUS BRANDS AND MODELS PHARMACEUTICAL REFRIGERATORS AT VACCINATION CENTERS BERAKAS HEALTH SCREENING.

1. BME7312, BRAND: PHILIP KIRSCH, MODEL: MED-125, SN: 1250101794 (Lampiran 2, item no. 2)
2. BME21272, BRAND: HAIER, MODEL: HYC-118A, SN: BE0691EBA00QEJ9H0001
3. BME5642, BRAND: SANYO, MODEL: MPR 513, SN: 41217572
4. BME3112, BRAND: SANYO, MODEL: MPR-411FR, SN: 90100023
5. BME24956, BRAND: HAIER, MODEL: HYC-290, SN: BE06M EE10 00EK9 A0003
6. BME24957, BRAND: HAIER, MODEL: HYC-290, SN: BE06M EE1T 00QEK 990019
7. BME24957, BRAND: HAIER, MODEL: HYC-290, SN: BE06M EE1T 00QEK 9A0044
8. BME5789, BRAND: PHILIP KIRSCH, MODEL: MED-340, SN: 3301000388
9. BME5790, BRAND: PHILIP KIRSCH, MODEL: MED-340, SN: 3301000391
10. BME5760, BRAND: PHILIP KIRSCH, MODEL: MED-520, SN: 5000510148 (Lampiran 2, item no. 3)
11. BME5757, BRAND: PHILIP KIRSCH, MODEL: MED-520, SN: 5000510138
12. BME5756, BRAND: PHILIP KIRSCH, MODEL: MED-520, SN: 5000510136 (Lampiran 2, item no. 4)

SCOPE OF WORK INCLUDING BUT NOT LIMITED TO;

- GENERAL SERVICE (INCLUDE CLEANING COOLING COIL, FAN BLOWER, DOOR SEALANT ETC)
- CLEAN AND CHECK THERMISTER
- RECALIBRATION TEMPERATURE SETTING
- TOP UP GAS (IF NECESSARY)
- INSPECT COMPRESSOR AND DOOR SEALANT.
- CLEAN BASE PLATE TO ENSURE STABILITY
- INSPECT CONTROL MONITORING DISPLAY UNIT.
- CHECK GREASE FAN MOTORS.
- PLEASE PROVIDE PROPER COPY OF PPM CHECK LIST.
- APPLY PPM STICKER/TAG AFTER SUCCESSFUL SERVICE.
- WORK ON PPM OF PHARMACEUTICAL REFRIGERATOR TO BE CARRIED OUT BY TRAINED ENGINEER(S).

LAMPIRAN - A2

RUJ: BME/A2K/(1194)/2023/HTD (B8)

ONE-OFF PREVENTIVE MAINTENANCE FOR TWELVE (12) UNITS OF VARIOUS BRANDS AND MODELS PHARMACEUTICAL REFRIGERATORS AT VACCINE CENTERS BERAKAS HEALTH SCREENING, MINISTRY OF HEALTH.

NO	SPAREPART	PART NUMBER	QUANTITY
1	SERVICE AND PPM		12
2	EVAPORATOR MOTOR ASSY		1
3	COMPRESSOR ASSY		1
4	THERMISTOR		1