Rujukan Kami: (46) MOH/HQ/P/IKLAN-SH/2024

BIL	Quotation Reference	Description	Advertisement Date	Closing Date (Not Later Than 02.00PM)	Quotation Fee	Requesting Department	Focal Person
4	(18) PPN/QTN/2024 (LAB)	SUPPLY AND DELIVERY OF LABORATORY CONSUMABLE ITEMS FOR DENTAL LABORATORY FOR DEPARTMENT OF DENTAL SERVICES, MINISTRY OF HEALTH PART (1)	26/03/2024	20/04/2024	\$5.00	JABATAN PERKHIDMATAN PERGIGIAN, PUSAT PERGIGIAN NEGERA, KEMENTERIAN KESIHATAN.	Hajah Siti Noriyah binti Haji Md Noor Senior Dental Technologists Department of Dental Services National Dental Centre Ministry of Health
	,	PLACE OF SUBMISSION: QUOTATION BOX (GROUND FLOOR) MINISTRY OF HEALTH COMMONWEALTH DRIVE BANDAR SERI BEGAWAN, BB 3910 NEGARA BRUNEI DARUSSALAM					

NO	ITEM DESCRIPTION	QUANTITY	BRAND	MANUFACTURER AND ORIGIN	UNIT PRICE	TOTAL PRICE
1.	Flasks	50 Units		Y		
	Specification for Flasks:					
	 For curing denture and appliances To hold plaster, model and acrylic inside Can fit in any clamp Preferably light weight In 2 parts – deep and shallow parts Packing: Pack of 1 					
2.	Clamp for flask Specification for Clamp for flask: To keep flasks held tightly	10 Units				
	during curing and prevent porosity during curing Able to fit normal type of flasks Screw on each rod to secure flasks in place With metal base and another metal piece (same size) Height: 20-26cm Packing: Pack of 1					

NO	ITEM DESCRIPTION	QUANTITY	BRAND	MANUFACTURER AND ORIGIN	UNIT PRICE	TOTAL PRICE
3.	Model Surveyor	3 Units	" # # # # # # # # # # # # # # # # # # #			
	Specification for Model Surveyor:		3 3	-		
	 Precision bearings at all joints 2 horizontal joints (double hinged arm) Twofold vertical adjustable working level Solid baseplate for 360° stability Detachable precision model table with secure model locking device Model table and tool holding fixture can be swivelled 360° Model holding plate can be swivelled with 0° lock Use of a second work table 					
	Comprehensive basic equipment Optimal sect benefit ratio					
	 Optimal cost benefit ratio Approximate dimension: 315mm (H) x 195mm (W) Approximate weight: 3.25kg 	-				

NO	ITEM DESCRIPTION	QUANTITY	BRAND	MANUFACTURER AND ORIGIN	UNIT PRICE	TOTAL PRICE
4.	Lab Putty Base	5 Buckets	П			
	Specification for Lab Putty Base:		(F			
	 Should mix well when used together with lab putty activator Adequate setting time – not more than 3 minutes Do not stick on denture or appliances Packing – bucket 7.5kg 					
5.	Lab Putty Activator Specification for Lab Putty Activator: • Should mix well when used together with lab putty base	50 Tubes				
	 Adequate setting time – not more than 3 minutes Do not stick on denture or appliances Longer shelf life (expiry date) – at least 18 months Packing – box of 40ml 		5			

NO	ITEM DESCRIPTION	QUANTITY	BRAND	MANUFACTURER AND ORIGIN	UNIT PRICE	TOTAL PRICE
6.	Thermoplastic Material for Essix Retainer	15 Boxes				
	 Specification for Thermoplastic Material for Essix Retainer: To be used with scheu vacuum machine Round shape – diameter 125mm x 1mm thick Non-mutagenic, non-toxic and non-irritating to mucous membrane Packing: Box of 100pcs 					

ASO/pg4/6

NO	ITEM DESCRIPTION	QUANTITY	BRAND	MANUFACTURER AND ORIGIN	UNIT PRICE	TOTAL PRICE
7.	Plaster of Paris Specification for Plaster of Paris: Must be dental plaster White colour For mounting, basing and flasking Setting time in the range of 12-20 mins Packing in the range of 20-25kg per bag Total amount required: 1400kg Shelf life of at least 18 months upon delivery Comply with ISO6873:2013 (Dentistry Gypsum Products) or with AS1652-1974 (Dental Laboratory Plaster)	70 Bags of at least 20kg/bag				

ASO/pg5/6

DELIVERY OF PERIOD: Not more than 60 days upon confirmation	(Yes/No) (If No, please specify)	
no supplier may withdraw his/her quotation	THS from the final date for the submission of the quotation and main with that period. The Government reserves the right to extend at such extension to the quotation validity period shall have the	

ASO/pg6/6