

BIL	Quotation Reference	Description	Advertisement Date	Closing Date (Not Later Than 02.00PM)	Quotation Fee	Requesting Department	Focal Person
5	(19) PPN/QTN/2024 (LAB)	SUPPLY AND DELIVERY OF LABORATORY CONSUMABLE ITEMS FOR DENTAL LABORATORY FOR DEPARTMENT OF DENTAL SERVICES, MINISTRY OF HEALTH PART (2) <u>PLACE OF SUBMISSION:</u> QUOTATION BOX (GROUND FLOOR) MINISTRY OF HEALTH COMMONWEALTH DRIVE BANDAR SERI BEGAWAN, BB 3910 NEGARA BRUNEI DARUSSALAM	26/03/2024	20/04/2024	\$5.00	JABATAN PERKHIDMATAN PERGIGIAN, PUSAT PERGIGIAN NEGERA, KEMENTERIAN KESIHATAN.	Hajah Siti Noriyah binti Haji Md Noor Senior Dental Technologists Department of Dental Services National Dental Centre Ministry of Health Tel : 2380433

QTN REF: (19) PPN/QTN/2024 (LAB)

**SUPPLY AND DELIVERY OF LABORATORY CONSUMABLE ITEMS FOR DENTAL LABORATORY, DEPARTMENT OF DENTAL SERVICES, MINISTRY OF HEALTH
PART (2)**

NO	ITEM DESCRIPTION	QUANTITY	BRAND	MANUFACTURER AND ORIGIN	UNIT PRICE	TOTAL PRICE
1.	Light-Cure Acrylic for Special tray Specification for Light-Cure Acrylic for Special tray: <ul style="list-style-type: none">• Allows easy fabrication of custom impression tray• Can be cured under halogen and UVA light• Minimum shrinkage and slump resistant• Non-sticky and easy to adapt and trim at the dough stage• Upper shape• Packaging: box of 50 boxes	50 Boxes				
2.	Stone Type 3 Specification for Stone Type 3: <ul style="list-style-type: none">• Must be a dental stone• Coloured – not white• Moderate setting time (10-12 mins)• Box of 20-25kg• Shelf life of at least 18 months upon delivery• Complies with ISO 6873:2013, Dentistry Gypsum Products	70 Boxes				

QTN REF: (19) PPN/QTN/2024 (LAB)

**SUPPLY AND DELIVERY OF LABORATORY CONSUMABLE ITEMS FOR DENTAL LABORATORY, DEPARTMENT OF DENTAL SERVICES, MINISTRY OF HEALTH
PART (2)**

NO	ITEM DESCRIPTION	QUANTITY	BRAND	MANUFACTURER AND ORIGIN	UNIT PRICE	TOTAL PRICE
3.	Resin Separator	20 Bottles				
	Specification for Resin Separator:					
	<ul style="list-style-type: none">• Used as a separating medium of plaster / resin and flask / plaster surfaces• Package: Bottle of 300g• Produces a durable surface film which does not peel off during trial closures• Gives easy separation of surfaces so that cured acrylic has smooth and shiny surface• Has long shelf life of at least 24 months upon delivery• Can also be used to coat inside of flasks before pouring for them to be emptied easily afterwards without distortions by prolonged hammerings					

QTN REF: (19) PPN/QTN/2024 (LAB)

**SUPPLY AND DELIVERY OF LABORATORY CONSUMABLE ITEMS FOR DENTAL LABORATORY, DEPARTMENT OF DENTAL SERVICES, MINISTRY OF HEALTH
PART (2)**

NO	ITEM DESCRIPTION	QUANTITY	BRAND	MANUFACTURER AND ORIGIN	UNIT PRICE	TOTAL PRICE
4.	Plastic Model for Orthodontic Study Models	10 Boxes				
	Specification for Plastic Model for Orthodontic Study Models: <ul style="list-style-type: none">• Consists of an upper base, a lower base and hinges• Bases should be shaped as a study model (angles trimmed)• Packing – box of 100 sets of upper and lower bases and 200pcs of hinges					
5.	Modelling Wax	100 Boxes				
	Specification for Modelling Wax: <ul style="list-style-type: none">• For dental use• Pink colour• Not a hard type of wax – moderate hardness• Does not distort in the mouth• Packing – box of 500g					

QTN REF: (19) PPN/QTN/2024 (LAB)

**SUPPLY AND DELIVERY OF LABORATORY CONSUMABLE ITEMS FOR DENTAL LABORATORY, DEPARTMENT OF DENTAL SERVICES, MINISTRY OF HEALTH
PART (2)**

NO	ITEM DESCRIPTION	QUANTITY	BRAND	MANUFACTURER AND ORIGIN	UNIT PRICE	TOTAL PRICE
6.	Duplicating Liquid Concentrate (for model):	20 Bottles of 2.7kg/bottle				
	Specification for Duplicating Liquid Concentrate (for model): <ul style="list-style-type: none">• Thermoreversible duplicating gel for prosthetics (duplicate models)• Shelf life of at least 24 months upon delivery• Total amount required: 54kg (To specific packaging size not more than 10kg per bottle)• Comply with ISO14356:2003, Dentistry Duplicating Material					

QTN REF: (19) PPN/QTN/2024 (LAB)

SUPPLY AND DELIVERY OF LABORATORY CONSUMABLE ITEMS FOR DENTAL LABORATORY, DEPARTMENT OF DENTAL SERVICES, MINISTRY OF HEALTH
PART (2)

<p>DELIVERY OF PERIOD:</p> <p>Not more than 60 days upon confirmation</p>	<p>(Yes/No) (If No, please specify)</p>	
<p>PRICE VALIDITY:</p> <p><i>The quotation shall remain valid for 6 MONTHS from the final date for the submission of the quotation and no supplier may withdraw his/her quotation with that period. The Government reserves the right to extend this period if deemed necessary provided that such extension to the quotation validity period shall have the written consent of the suppliers.</i></p>		

ASO/pg5/5

QTN REF: (19) PPN/QTN/2024 (LAB)

**SUPPLY AND DELIVERY OF LABORATORY CONSUMABLE ITEMS FOR DENTAL LABORATORY,
DEPARTMENT OF DENTAL SERVICES, MINISTRY OF HEALTH PART (2)**

NO	ITEM DESCRIPTION	USER REQUIREMENT	YES	NO	PLEASE SPECIFY THE SPECIFICATION / SCOPE OFFERED FAILURE TO FILL THE OFFER DETAIL SHALL BE IDENTIFIED AS NON-COMPLIANCE TO THE OFFER
1.	LIGHT-CURE ACRYLIC FOR SPEACIAL TRAY	Allows easy fabrication of custom impression tray			
		Can be cured under halogen and UVA light			
		Minimum shrinkage and slump resistant			
		Non-sticky and easy to adapt and trim at the dough stage			
		Upper shape			
		Packaging: box of 50 boxes			
		50 Boxes			
2.	STONE TYPE 3	Must be a dental stone			
		Coloured – not white			
		Moderate setting time (10-12 mins)			
		Box of 20-25kg			
		Shelf life of at least 18 months upon delivery			
		Complies with ISO 6873:2013, Dentistry Gypsum Products			
		70 Boxes			

QTN REF: (19) PPN/QTN/2024 (LAB)

**SUPPLY AND DELIVERY OF LABORATORY CONSUMABLE ITEMS FOR DENTAL LABORATORY,
DEPARTMENT OF DENTAL SERVICES, MINISTRY OF HEALTH PART (2)**

NO	ITEM DESCRIPTION	USER REQUIREMENT	YES	NO	PLEASE SPECIFY THE SPECIFICATION / SCOPE OFFERED FAILURE TO FILL THE OFFER DETAIL SHALL BE IDENTIFIED AS NON-COMPLIANCE TO THE OFFER
3.	RESIN SEPARATOR	Used as a separating medium of plaster / resin and flask / plaster surfaces			
		Package: Bottle of 300g			
		Produces a durable surface film which does not peel off during trial closures			
		Gives easy separation of surfaces so that cured acrylic has smooth and shiny surface			
		Has long shelf life of at least 24 months upon delivery			
		Can also be used to coat inside of flasks before pouring for them to be emptied easily afterwards without distortions by prolonged hammerings			
		20 Bottles			
4.	PLASTIC MODEL FOR ORTHODONTIC STUDY MODELS	Consists of an upper base, a lower base and hinges			
		Bases should be shaped as a study model (angles trimmed)			
		Packing – box of 100 sets of upper and lower bases and 200pcs of hinges			
		10 Boxes			

QTN REF: (19) PPN/QTN/2024 (LAB)

**SUPPLY AND DELIVERY OF LABORATORY CONSUMABLE ITEMS FOR DENTAL LABORATORY,
DEPARTMENT OF DENTAL SERVICES, MINISTRY OF HEALTH PART (2)**

NO	ITEM DESCRIPTION	USER REQUIREMENT	YES	NO	PLEASE SPECIFY THE SPECIFICATION / SCOPE OFFERED FAILURE TO FILL THE OFFER DETAIL SHALL BE IDENTIFIED AS NON-COMPLIANCE TO THE OFFER
5.	MODELLING WAX	For dental use			
		Pink colour			
		Not a hard type of wax – moderate hardness			
		Does not distort in the mouth			
		Packing – box of 500g			
		100 Boxes			
6.	DUPLICATING LIQUID CONCENTRATE (FOR MODEL)	Thermoreversible duplicating gel for prosthetics (duplicate models)			
		Shelf life of at least 24 months upon delivery			
		Total amount required: 54kg (to specific packaging size not more than 10kg per bottle)			
		Comply with ISO14356:2003, Dentistry Duplicating Material			
		20 Bottles of 2.7kg/bottle			