

LAMPIRAN 1

BIL	Quotation Reference	Description	Advertisement Date	Closing Date (Not Later Than 02.00PM)	Quotation Fee	Requesting Department	Focal Person
1	BME/A2K/(1253)/2024/HTD(A20)	<p>ONE-OFF PREVENTIVE MAINTENANCE FOR ONE UNIT (1) OF GAS CHROMATOGRAPHY MASS SPECTROMETRY (GCMS) SYSTEM, BRAND: SHIMADZU MODEL: QP2010 PLUS, S/N: C70504501120, BME6686 AT FOOD CHEMISTRY SECTION, DEPARTMENT OF SCIENTIFIC SERVICES</p> <p><u>PLACE OF SUBMISSION:</u> QUOTATION BOX (GROUND FLOOR) MINISTRY OF HEALTH COMMONWEALTH DRIVE BANDAR SERI BEGAWAN, BB 3910 NEGARA BRUNEI DARUSSALAM</p>	23/05/2024	08/06/2024	\$5.00	JABATAN TEKNOLOGI PENJAGAAN KESIHATAN, KEMENTERIAN KESIHATAN.	<p>Haji Zaenul 'Arifien bin Dato Paduka Haji Mohd Rosli Jurutera Perubatan Hayat Jabatan Teknologi Penjagaan Kesihatan</p> <p>Tel : 2382424 Ext : 6122/6232</p>

LAMPIRAN - A1

RUJ: BME/A2K/(1253)/2024/HTD (A20)

9/5/2024

ONE-OFF PREVENTIVE MAINTENANCE FOR ONE UNIT (1) OF GAS CHROMATOGRAPHY MASS SPECTROMETRY (GCMS) SYSTEM, BRAND: SHIMADZU MODEL: QP2010 PLUS, S/N: C70504501120, BME6686 AT FOOD CHEMISTRY SECTION, DEPARTMENT OF SCIENTIFIC SERVICES.

SCOPE OF WORK

One-off Preventive Maintenance for one unit (1) of Shimadzu GCMS-QP2010 PLUS, S/N: C70504501120, BME6686 at Food Chemistry Section, Department of Scientific Services, inclusive of but not limited to:

- Scope of work to follow preventive maintenance manual/recommendation specific for this model, which include the following where necessary:

- (i) inspection
- (ii) cleaning
- (iii) installation of PM kits, consumables and spare parts listed in spareparts list (Lampiran A2)
- (iv) alignment
- (v) calibration
- (vi) any other related preventive maintenance works required.

- Functional check post-maintenance.

- Apply preventive maintenance sticker/tag after successful service.

- Disposal of any used spare parts/pm kits, unless otherwise instructed.

- Work on the system is to be carried out by engineer(s) authorised in providing maintenance to the equipment.

- Vendor(s) to submit maintenance checklist/details/scope specific for this equipment along with other necessary documents to join this tender.

- Upon approval, the earliest the preventive maintenance work can be carried out on the equipment is one month before the next PM due date, and not anytime earlier.

LAMPIRAN - A2

RUJ: BME/A2K/(1253)/2024/HTD (A20)

ONE-OFF PREVENTIVE MAINTENANCE FOR ONE UNIT (1) OF GAS CHROMATOGRAPHY MASS SPECTROMETRY (GCMS) SYSTEM, BRAND: SHIMADZU MODEL: QP2010 PLUS, S/N: C70504501120, BME6686 AT FOOD CHEMISTRY SECTION, DEPARTMENT OF SCIENTIFIC SERVICES.

NO	SPAREPART	PART NUMBER	QUANTITY
1	GCMS PM Kit	225-02091-41	1
2	Rotary pump oil ultragrade 15, 1L/btl	017-30163-11	1
3	Ceramic Bush for Insulator	225-01068	4
4	Ceramic Insulator	200-44394	10
5	Split insert with wool (5's/pkt)	221-75193	1
6	Grindstine, Hp13S, M2 (12's/pkt)	670-18701-02	1
7	Shimadzu 1 Triple Helium specific Filter - 1/8"SS	226-50712-00	1
8	SH-RTx-5MS Cap. Column ID 0.25mm, 0.25um, 30m	221-75855-30	1