

BIL	Quotation Reference	Description	Advertisement Date	Closing Date (Not Later Than 2.00PM)	Quotation Fee	Requesting Department	Focal Person
5	BME/A2K/(995)/2023/HTD (A20)	<p>ONE OFF PREVENTIVE MAINTENANCE FOR FOUR (4) UNITS OF HYDROGEN GAS GENERATOR, BRAND: PARKER, MODEL: GCGS 7890-220, AT DEPARTMENT OF SCIENTIFIC SERVICES</p> <p>PLACE OF SUBMISSION: QUOTATION BOX (GROUND FLOOR) MINISTRY OF HEALTH COMMONWEALTH DRIVE BANDAR SERI BEGAWAN, BB 3910 NEGARA BRUNEI DARUSSALAM</p>	06/02/2023	25/02/2023	\$5.00	JABATAN TEKNOLOGI PENJAGAAN KESIHATAN, KEMENTERIAN KESIHATAN.	ZAENUL 'ARIFIEN BIN DATO PADUKA HAJI MOHD ROSLI JURUTERA PERUBATAN HAYAT TEL: 2382424 ext. 6122 / 6232

LAMPIRAN - A1

RUJ: BME/A2K/(995)/2023/HTD (A20)

2/2/2023

ONE OFF PREVENTIVE MAINTENANCE FOR FOUR (4) UNITS OF HYDROGEN GAS GENERATOR, BRAND: PARKER, MODEL: GCGS 7890-220, AT DEPARTMENT OF SCIENTIFIC SERVICES.

SCOPE OF WORK

One off preventive maintenance services for four (4) units of hydrogen gas generators, Brand: Parker, Model: GCGS 7890-220, at Department of Scientific Services.

The details of each generator are as follows:

1. S/N: GCGS7890220011, at Food Chemistry Section.
2. S/N: GCGS7890220013, at Food Chemistry Section
3. S/N: GCGS7890220010, at Illicit Drugs Section.
4. S/N: GCGS7890220018, at Toxicology Section.

- These four units each require the 12 month maintenance kit, where each kit must consist of:

- a. Hydrogen Filter (DI Resin Cartridge) (P/N: B02-0503)
- b. Membrane Filter (Zero Air Outlet Filter) (P/N: GS050-05-95)
- c. Coalescing Filter (Zero Air Intlet Filter) (P/N: 050-05-BX)

- Scope of work on the hydrogen gas generators to follow Parker preventive maintenance manual/recommendation specific for model GCGS7890220, which include the following where necessary:

- (i) inspection
- (ii) cleaning
- (iii) installation of PM kits, consumables and spare parts listed in spareparts list (Lampiran A2)
- (iv) alignment
- (v) calibration
- (vi) any other related preventive maintenance works required

- Functional check post-maintenance.

- Apply preventive maintenance sticker/tag after successful service.

- Disposal of any used spare parts/pm kits, unless otherwise instructed.

- Work on the hydrogen gas generators to be carried out by authorised engineer(s) authorised in providing maintenance to Parker gas generators.

LAMPIRAN - A2

RUJ: BME/A2K/(995)/2023/HTD (A20)

ONE OFF PREVENTIVE MAINTENANCE FOR FOUR (4) UNITS OF HYDROGEN GAS GENERATOR, BRAND: PARKER, MODEL: GCGS 7890-220, AT DEPARTMENT OF SCIENTIFIC SERVICES.

NO	SPAREPART	PART NUMBER	QUANTITY
1	12 Month Maintenance Kit		4